



Data-Driven Healthcare at current health

Annabel Kramer

Data Scientist

currenthealth.com





@current_health



snap40 is now current health!

Outline


- Introduce our device & platform
- Device data
- Challenges & considerations when working with health data
- Published paper: Recognising Cardiac Abnormalities in Wearable Device PPG Waveforms using Deep Learning






Temperature

- Skin & ambient temperature sensors
- Detecting hypo-/hyperthermia





https://github.com/nlathia/pydata_2016/blob/master/Talk/Presentation.pdf

- Triaxial accelerometer and gyroscope
- Useful for contextualising patient vital signs + behaviour
- Many more applications:
 - Fall detection
 - Consciousness/sleep
 - Shiver/cough detection



Respiration Rate
+/- 1 rpm



Pulse Rate
+/- 3bpm



Oxygen Saturation
+/- 2%



Temperature
+/- 0.1°C



**Movement and
Posture**

What's the heart rate of a banana?




<https://www.macobserver.com/news/check-banana-heart-rate-apple-watch/>

Considerations & Challenges

- Provide accurate, meaningful information
 - Reduce alarm fatigue: less false positives
 - Actionable information
- Confounding events
 - What else could have affected what I'm looking at?
- Validation
 - Medical gold standard
 - Population bias & test set distribution
 - Regulatory requirements

Validation: FDA




Recognising Cardiac Abnormalities in Wearable Device PPG with Deep Learning


Identify whether a PPG waveform contains any cardiac abnormalities

<https://arxiv.org/pdf/1807.04077.pdf>


PPG Data




<https://www.mdpi.com/2079-9292/3/2/282/htm>




PPG Data




(a)



(b)



(c)




(d)


(a) shows a normal PPG waveform, (b)-(d) contain abnormalities

Method


- Train an LSTM autoencoder to reconstruct PPG waveforms
- Measure agreement between original and reconstructed waveform
- Low agreement → abnormality
- Evaluate approach on cross-section of hospital and at-home patients

—●— Input PPG signal -x- Auto-encoded PPG signal







(a)



(b)



(c)



(d)

- 60% true positives (anomaly detected in PPG, and anomaly present in ECG)
- 23% false positives (anomaly detected in PPG but *not* present in ECG)
- Promising results
- Can we cluster together similar anomalies?



Thank you!

Questions?

currenthealth.com



[@current_health](https://twitter.com/current_health)



Want to get in touch?

Data Science: annabel@snap40.com

Sales: rachel@snap40.com

Hiring: people@snap40.com

Check out our website snap40.com

New website to launch soon!

currenthealth.com



@current_health